

## imaet PRO Version 11.3 Release Notes

Release Date: 1/21/24

### New Software Features & Improvements

- Category Button glitch on Biofield Panel fixed
- New Items and Item Descriptions added to the Allergen Panel

### Box and Accessories

- None

### Minimum System Requirements

- Windows 10 or later
- 8 Gigabytes RAM
- i5 Processor
- 512 Gigabytes of Internal Memory

### Known Issues & Limitations

- The Biofield Panel is intended for advanced users. Newer users should utilize the Allergen Panel instead until they have completed QBA training or after a minimum of 6 months of use.
- Geopathic Panel has a visual glitch involving the Interface button continuing to blink after the feedback has completed. Clicking Stop Feedback will stop the blinking. This does not affect functionality.
- New Chart function on Client Profile not yet functional.
- Panel Info page is unavailable on the Pet Panel.
- Anti-virus or anti-malware software running during use can cause malfunction.
- Pending Windows Updates and driver updates can cause malfunction.
- Opening the program while it is already running, creating multiple instances of the program, causes malfunction.

See our Recommended [Device Best Practices](#) for optimal performance